

pharmbio robotlab


Dan Rosén, 2023-08-30

High-content cell profiling


robotlab cellpainter

robotlab imager


① Cell seeding + treatment




② Cell Painting




③ Image acquisition



④ Feature extraction



⑤ Data analysis



Cell painting protocol

Each plate goes through:


- 5 washes (70s-135s)
- 4 dispenses (20s-30s)
- 4 incubation (20 min)

Lid on when incubating

Lid off when using washer and dispenser

Filter (excitation/emission)	Exposure ¹ (ms)	Dye	Organelle or cellular component	Channel name, in CellProfiler
DAPI (387/447 nm)	100	Hoechst 33342	Nucleus	DNA
GFP (472/520 nm) ²	100	concanavalin A (con A) AlexaFluor488 conjugate	Endoplasmic reticulum ²	ER
Cy3 (531/593 nm)	200	SYTO 14 green fluorescent nucleic acid stain	Nucleoli, cytoplasmic RNA ³	RNA
TexasRed (562/624 nm) ⁴	100	Alexa Fluor 594 phalloidin conjugate, wheat germ agglutinin (WGA) AlexaFluor594 conjugate	F-actin cytoskeleton, Golgi, plasma membrane	AGP
Cy5 (628/692 nm)	400	MitoTracker Deep Red	Mitochondria	Mito


Cell Painting, an image-based assay for morphological profiling (Bray et al, 2016)



Automated cell painting



37C incubator
(Liconic)



Robot arm (UR)



Biotek dispenser



Biotek washer

Components

github.com/pharmbio/robotlab

- Machine communication [labrobots/](#)
 - HTTP server (flask) wrapping lower-level communication (COM-ports, subprocesses, ...)
 - Type-safe python wrapper
 - Developers can access using curl
- Robotarm moves [cellpainter/{moves_gui.py,movelist/*jsonl}](#)
 - Custom UI
 - Interpolate between hotel positions
 - Moves cut up into smaller pieces using code
- Scheduling [cellpainter/{commands,optimize,protocol}.py](#)
 - DSL for machine execution and time checkpointing.
 - SMT solver Z3 to solve waiting delays.
- Front end web browser UI [cellpainter/gui/](#)
 - Generated server-side from python
 - CLI also present but never used

Fridge communication



```
dan@devserver:~$ curl 10.10.0.97:5050
{
  "http://10.10.0.97:5050/echo": "Echo()",
  "http://10.10.0.97:5050/git": "Git()",
  "http://10.10.0.97:5050/fridge":
    "Fridge(id='STX', host='localhost', port=3333, mode='execute',
      current_climate=Cell(value={'temp': 4.5, 'humid': 83.6, 'co2': 0.0, 'n2': 240.23})),
      fridge_db='fridge.db')",
  "http://10.10.0.97:5050/barcode": "BarcodeReader(com_port='COM3', current_barcode=Cell(value='P102813'))",
  "http://10.10.0.97:5050/io.sql": "Download the IO database as sqlite dump",
  "http://10.10.0.97:5050/io.db": "Download the IO database in binary sqlite",
  "http://10.10.0.97:5050/tail": "Show last 10 lines from the IO database",
  "http://10.10.0.97:5050/tail/<N>": "Show last N lines from the IO database"
}
```

```
dan@devserver:~$ curl 10.10.0.97:5050/tail
280319 2023-08-29 13:31:03.793 fridge 48879 eject(plate='P102813', project='SSS-val')
280320 2023-08-29 13:31:03.808 fridge 48879 stx.write(b'STX2ServiceMovePlate(STX,2,3,2,1,1,STX,1,0,0,1,1)\r')
280321 2023-08-29 13:31:35.048 fridge 48879 stx.read() = b'1\r\n'
280322 2023-08-29 13:31:35.115 fridge 48879 31312.0ms eject(plate='P102813', project='SSS-val')
280323 2023-08-29 13:31:35.124 fridge 48879 return ('3x2', {'plate': 'P102813', 'project': 'SSS-val'})
```

Fridge communication

```
@dataclass(frozen=True)
class Fridge(STX):
    fridge_db: str = 'fridge.db'

    def eject(self, plate: str, project: str) -> tuple[str, FridgeSlot]:
        """
        Ejects a plate given a plate and project.

        Returns info about the plate and the old location,
        or raises an error if it was not possible to complete the action.
        """
        with self.exclusive():
            with self._get_db() as db:
                loc_slot = db.get_by_plate_project(plate, project)
                if not loc_slot:
                    raise ValueError(f'No slot with {plate=} and {project=}')
                loc, slot = loc_slot
            return self._eject(loc, slot)
```

```
dan@devserver:~$ curl 10.10.0.97:5050/fridge
{
  "http://10.10.0.97:5050/fridge/eject": [
    "eject(plate: 'str', project: 'str') -> 'tuple[str, FridgeSlot]'",
    "",
    "Ejects a plate given a plate and project.",
    "",
    "Returns info about the plate and the old location,",
    "or raises an error if it was not possible to complete the action."
  ],
  ...
}
```

HTTP endpoints generated from the class

```
curl .../fridge/eject/P001337/protac
```

```
curl -X POST .../fridge/eject -d '{
  "plate": "P1337", "project": "protac"}
```

Fridge communication

```
@dataclass(frozen=True)
class Fridge(STX):
    fridge_db: str = 'fridge.db'

    def eject(self, plate: str, project: str) -> tuple[str, FridgeSlot]:
        '''
        Ejects a plate given a plate and project.

        Returns info about the plate and the old location,
        or raises an error if it was not possible to complete the action.
        '''
        with self.exclusive():
            with self._get_db() as db:
                loc_slot = db.get_by_plate_project(plate, project)
                if not loc_slot:
                    raise ValueError(f'No slot with {plate=} and {project=}')
                loc, slot = loc_slot
            return self._eject(loc, slot)
```

```
class Runtime:
    def __post_init__(self):
        gbg = WindowsGBG.remote()
        self.fridge = gbg.fridge
        self.barcode_reader = gbg.barcode

    def execute(cmd: Command, runtime: Runtime):
        match cmd:
            case FridgeEject():
                runtime.fridge.eject(
                    plate=cmd.plate, project=cmd.project)
```

Type-safe communication:

Type-checker thinks `runtime.fridge` is a `Fridge`,
but it is a proxy object that does RPC over HTTP.

GUI for creating robot arm moves

Index x +
http://localhost:5000/?program=wash§ion=put

B21 disp h11 incu incu_A21 lid_h19 neu out21 r21 **wash** wash21 wash_to_disp wash_to_r21

wash put

run from here	MoveJoint(94.6, -111.3, 93.2, 18.0, 94.6, -0.1)	update			go	h neu
run from here	MoveLin(198.7, -574.7, 809.1, 0.0, 0.0, 90.0)	update	0.1	-5.8	0.0	go h21 neu
run from here	MoveLin(198.7, -574.6, 790.4, 0.0, 0.0, 90.0)	update	0.1	-5.7	-18.7	go h21 pick
run from here	GripperMove(255)	update			gripper close	go
run from here	MoveLin(198.7, -574.6, 809.1, 0.0, 0.0, 90.0)	update	0.1	-5.7	0.0	go h21 neu
run from here	MoveLin(198.7, -421.7, 809.1, 0.0, 0.0, 90.0)	update	0.1	147.2	0.0	go h neu
run from here	MoveLin(-2.0, -422.2, 808.2, 0.0, 0.0, 55.1)	update	-200.6	146.7	-0.9	go h safe
run from here	MoveJoint(69.7, -114.5, 95.1, 19.4, 104.6, -0.1)	update				go h safe
run from here	MoveJoint(-96.6, -126.2, 147.6, -22.9, -5.4, 1.5)	update				go wash above
run from here	MoveLin(-269.0, 195.2, 347.7, 0.0, 0.0, -1.2, slow=True)	update			1008	go wash above
run from here	MoveLin(-269.0, 195.2, 148.4, 0.0, 0.0, -1.3, slow=True)	update			1113	go wash neu
run from here	MoveLin(-269.0, 195.2, 132.9, 0.0, 0.0, -1.3, slow=True)	update			1122	go wash drop
run from here	GripperMove(89)	update			gripper open	go
run from here	MoveLin(-269.0, 195.2, 148.4, 0.0, 0.0, -1.3, slow=True)	update			1113	go wash neu
run from here	MoveLin(-269.0, 195.2, 347.7, 0.0, 0.0, -1.2, slow=True)	update			1008	go wash above
run from here	MoveJoint(-96.6, -126.2, 147.6, -22.9, -5.4, 1.5)	update				go wash above
run from here	MoveJoint(69.7, -114.5, 95.1, 19.4, 104.6, -0.1)	update				go h safe
run from here	MoveLin(-2.0, -422.2, 808.2, 0.0, 0.0, 55.1)	update	-200.6	146.7	-0.9	go h safe
run from here	MoveLin(198.7, -421.7, 809.1, 0.0, 0.0, 90)	update	0.1	147.2	0.0	go h neu

run program freedrive **stop robot** gripper open gripper close grip test



roll -> 0° (level roll) pitch -> 0° (face horizontally) pitch -> -90° (face the floor) yaw -> 0° (towards washer and dispenser) yaw -> 90° (towards hotels and incu)

```
{'xyz': [198.6, -568.9, 809.1],  
'rpy': [0.0, 0.0, 90],  
'joints': [93.06, -98.56, 84.14, 14.42, 93.06, 0.0],  
'pos': [88],  
'tick': [8],  
'server_age': 51}
```


set speed to 20
set speed to 40
set speed to 60
set speed to 80
set speed to 100

- Written for the Universal Robots robot, later ported to the PreciseFlex robot


Designing and sequencing robot arm moves



- We use a neutral position by the hotel as a reference point for all moves
- When sequencing two moves we can remove redundant returns to the neutral position
- We generate moves to all hotel positions by offsetting a chosen subset of the moves with the hotel separation distance



Expressing the protocol in code and scheduling it



The protocol is expressed in a domain-specific language (embedded in python):

```
cmd ::=      RunRobotarm(program_name)
           |  RunDispenser(protocol_name)
           |  Checkpoint(checkpoint_name)
           |  WaitForCheckpoint(checkpoint_name, plus_seconds)
           |  Fork(cmd)
           |  Sequence(cmd, cmd)
```

The program can be executed in three ways:

	1. Live	2. Simulation	3. Symbolically
Time:	Wall time	Simulated	Symbolic
Machines:	Execute	Advance time	Advance time
Purpose:	Side effects	Testing	Scheduling waits


Solve and
optimize using
off-the shelf
SMT-solver (Z3)

Example: cell painting a batch of two plates

Example: cell painting a batch of two plates

robot arm
dispenser
delay

batch [3:59.0]
incu [2:00.0, 2:00.0]



Protocol expressed in our Domain-Specific Language (DSL) embedded in Python.


```
cmd ::= RunRobotarm(program_name)
      | RunDispenser(protocol_name)
      | Checkpoint(checkpoint_name)
      | WaitForCheckpoint(checkpoint_name, plus_seconds)
      | Fork(cmd)
      | Sequence(cmd, cmd)
```

Scheduling by translation to SMT (fragment: Quantifier-Free Linear Arithmetic).

Example non-linearity


Case $|PFA| < |Mito|$:

batch [3:59.0]
incu [2:00.0, 2:00.0]



Case $|PFA| > |Mito|$:

batch [4:14.0]
incu [2:00.0, 2:00.0]



The delays have to be allocated differently depending on which case we are in. Formally, this means that the optimization problem is not linear.

Expressing complex protocols

Four plates in the MitoTracker cycle of the cell painting protocol.

